


ALPHACAM™ 1CMOS HD Camera [HD Vision-1C]

- 1/3" 1CMOS Super HD
- 1080P full HD Image output, life like colour reproduction
- Image freeze & 2.5 x Zoom function
- Highest resolution of 1000 lines for better depth & sharpness
- Camera head with four operatory functions
- Multiple output offer flexibility for connecting to peripheral Devices HDMIx1, CVBSx2, S-video, RGB
- Full HD recording in 1080P, Record Still picture & Live video into USB disk
- Anti-moire function for flexible endoscope



DEVICE

Part No.	Description
HD Vision-1C	Alphacam™ 1CMOS Full HD Camera with Full HD Video Recording with All Its Standard Accessories
HDV-1C-CH	1CMOS camera head (HD Vision-1C)
HDV-1C-CCU	1CMOS control unit (CCU) (HD Vision-1C)
ZCO-3	Zoom coupler 18-36mm

SUPPLIED ACCESSORIES

Part No.	Description
BNC-BNC	Monitor Cable
HDMI-DVI-D	HDMI to DVI-D cable
PD-16	Pen Drive 16GB
MC-03	Mains cord


TECHNICAL DATA

Camera Head	1/3' 1CMOS high brightness chip
Resolution	1920 x 1080P
High definition output (HD)	1080P Full HD
Minimum illumination	5 lux
Signal to Noise ratio	60db
Shutter	1/60s-1/10000s(automatic)
White Balance	Manual
Video output	HDMI, RGB, CVBSx2, S-VIDEO
Camera Head with control button	Yes with IPX8 waterproof grade
Camera Cable Length	3 meters
Power Supply	100-240 AC
Frequency	50 to 60 Hz
Maximum Power	50VA
Classification	Protection Class I, Type BF
Dimensions (WxHxD)	300 x 125 x 355 mm
Weight	4.5 Kgs.
Recording	Videos & Images (mpeg-4)
Video and image (formats)	MPEG & JPEG

Monitor Cable
BNC-BNC



*Composite
BNC-RCA



HDMI to DVI-D cable
HDMI-DVI-D



* HDMI to HDMI cable
HDMI-HDMI



Pen Drive 16GB
PD-16



Video Coupler Zoom [OMS],
Focal length : 18-36mm.
ZCO-3



Mains cord
MC-03



*Endoscope Coupler
Fixed Focus f=25mm
VCO-25

