


**ALPHFLATOR® CO<sub>2</sub> [AF-25L]**

- High safety
- Electronic controlled gas refilling technology
- Allows easy setup and monitoring of all parameters during insufflation process
- 0.1 to 25 l/min Flow
- Micro controller based
- Inbuilt heating system
- Pressure relief system



**DEVICE**

Part No.	Description
AF-25L	Alphflator® Insufflator for Laparoscopy, Co2 Gas 25 lpm LCD Display With Hose Pipe

**SUPPLIED ACCESSORIES**

Part No.	Description
ST-LL 3	Tube set, Insufflation reusable with luer lock
S-28/ S-16	Spanner
MC-03	Mains cord
HPT-04G	High Pressure tubing Co2 with guage

**TECHNICAL DATA**

Insufflation medium	Medical Carbon dioxide gas (Co2)
Display	LCD
Gas flow Maximum	25lpm
Pressure range	5-30mmHg
Power supply	100-240 V ac
Frequency	50-60Hz
Maximum Power	100VA
Classification	Protection Class 1, type Bf
Dimension (W x H X D)	300 x 138 x 320 mm
Weight	12kg

Tube set, Insufflation, reusable with luer lock  
ST-LL 3



Spanner  
S-28/ S-16



Mains cord  
MC-03



High Pressure tubing Co2 with guage  
HPT-04G



\* High Pressure Tubing for central gas supply- HPT-CS3



\* Filter single use only  
CF-01



\* High Pressure tubing Co2 pin index  
HPT-PI4G



Further Accessories available on request