


ALPHFLATOR® CO2 [AF-40LC]

- High safety
- Electronic controlled gas refilling technology
- Allows easy setup and monitoring of all parameters during insufflation process
- 0.1 to 40 l/min Flow
- Micro controller based
- Inbuilt heating system
- Pressure relief system
- Cylinder pressure display



DEVICE

Part No.	Description
AF-40LC	Alphflator® Insufflator for Laparoscopy, Co2 Gas 40 lpm LCD Display With Hose Pipe

SUPPLIED ACCESSORIES

Part no.	Description
HPT-45	High pressure tubing Co2
ST-LL 3	Tube set, insufflation reusable with luer lock
S-28/ S-16	Spanner
MC-03	Mains cord


TECHNICAL DATA

Insufflation medium	Medical Carbon dioxide gas (Co2)
Display	LCD
Gas flow Maximum	40 lpm
Pressure range	5-30mmHg
Power supply	100-240 V ac
Frequency	50-60Hz
Maximum Power	100 VA
Classification	Protection Class 1, type Bf
Dimension (W x H X D)	300 x 138 x 320 mm
Weight	12kg

Tube set, Insufflation , reusable with luer lock- ST-LL 3



* High Pressure Tubing for central gas supply- HPT-CS3



Spanner S-28/S-16



* Filter single use only CF-01



High pressure tubing co2 HPT-45



* High Pressure tubing Co2 pin index HPT-PI4



Mains cord MC-03



Further Accessories available on request