

ALPHKUT® ACCUSEAL PLUS [AA-150CB]

- Permanent and safe sealing up to 7mm.
- Seal can bear high blood pressure
- Resection can be performed with saline.
- Irrigation fluid remains clean for longer period.
- Saline virtually eliminates the risk of TUR Syndrome.
- Feather touch keyboard
- Bipolar cut function for sealing and cutting simultaneously
- LCD Graphics display
- Short operating time
- Isolated RF outputs
- Two modes for sealing pulsating and continuous.
- Lesser charring and drag.
- Automatic output adjustment dependent on tissue Impedance.



DEVICE

Part No.	Description
AA-150CB	Alphkut® Accuseal Plus [Complete Bipolar] [Electrosurgical Unit] With All Its Standard Accessories 150W (Bipolar,Bi-cut, Bi-Tur & Sealer)

SUPPLIED ACCESSORIES

Part no.	Description
FSC-03-AAP	Tri Pedal foot switch
NLPG-05-30-31-SA	Lap-Bi-Clamp Grasper
NLCP-05-30-31-SA	Lap Bi Clamp Cutter
MC-03	Mains Cord
BCS-04-03-VMF	Bipolar cable x2

TECHNICAL DATA

Bipolar			
Micro	80 Watt	100 OHM	3.0 CF
Macro	99 Watt	100 OHM	3.5 CF
Fast Seal L-Bip(L1-L4)	180 Watt	100 OHM	1.5 CF
Bi-Seal			
Pulsating L1-L4	120 Watt	100 OHM	1.5 CF
Continuous S1-S4	150 Watt	100 OHM	1.5 CF
Bi-Cut			
Cut (10-90)	200 Watt	100 OHM	1.5 CF
Bi-tur Cut			
Pure Cut	300 Watt	100 OHM	1.5 CF
Blend cut	250Watt	100 OHM	1.7 CF
Bitur Coag			
Spray	120Watt	100 OHM	4.5 CF
Fulguration	150Watt	100 OHM	4.0 CF
Dessication	120Watt	100 OHM	3.5 CF
Power supply	170 - 260 VAC		
Frequency	50 Hz		
Maximum Power	400W		
Classification	Protection Class 1, type Bf		
Dimension (W x H X D)	300 x 130 x 365 mm		
Weight	8kg		

Tri Pedal Footswitch
FSC-03-AAP



Lap-Bi Clamp Grasper
NLPG-05-30-31-SA



Lap-Bi Clamp Cutter
NLCP-05-30-31-SA



Mains cord
MC-03



Bipolar cable x2
BCS-04-03-VMF



*Open Bipolar Clamp
OBC-08



* Open Bipolar Clamp Cable
BCSO-04-03-2MF



*Resectoscope Monopolar/ Bipolar
MBPRC-4/302



Further Accessories available on request