

MEDICAL GRADE MONITOR

- TFT displays with wide screen format
- Back light stabilization system
- Built in multiple display mode
- True color restore
- Wide viewing angle
- Hygenic design with full glass front



Mains cord
MC-03



DEVICE

Part No.	Description
LM-24	Medical Grade Lemke Monitor 24"
LM-26	Medical Grade Lemke Monitor 26"
LM-27	Medical Grade Lemke Monitor 27"
MDSC-2226	Medical Grade Barco Monitor 26"
MDSC-2227	Medical Grade Barco Monitor 27"
P2415Q	Dell Monitor 24" [4K Monitor]
U2412M	Dell Monitor 24"

TECHNICAL DATA	LEMKE	BARCO
Model	LM27	MDSC 2226
Display size	27"	26"
Support color	1.07 billion	1 billion
Resolution	1920 x 1080(Full HD)	2MP(1920 x 1080)
Pixel Pitch	0.3114mm	0.3 mm
Brightness	700cd/m2 Typ.	450cd/m2 Typ.
Contrast	1000:1 Typ	1400:1 Typ
Angle of view	178° Typ	178° Typ
Response time	5ms. Typ	5ms. Typ
Input signal interface	DVI-D:24pin/3G-SDI: BNCYPbPr: BNC*3/RGBS:BNC*4 CVBS:BNC/ S-VIDEO: DIN4, VGA: D-Sub 15pin/ Fiber: FC(Optional)	DVI-I Single link (with HDMI support) RGBS/YPbPr (4xBNC) S-video (4-pin Mini DIN) Composite video (1xBNC) 3G-SDI (1xBNC)
Output signal interface	DVI-D:24pin/3G-SDI:BNC, CVBS:BNC, S-VIDEO:DIN4, YPbPr:BNC*3, RGBS:BNC*4	DVI-D x1 S-video (4-pin Mini DIN) Composite video (1xBNC) 3G-SDI (1xBNC)
Size	645.25(W)x392(H)x50(D)MM	633(W)x411(H)x60MM
Weight	9.0kgs.	9.4kgs
Weight with stand	12.5kg	12.5 kg
Standard Mounting	VESA 100x100mm	VESA (100mm)
Power supply	DC 24V	DC 24V
Power Consumption	55W max	120W max
Power Saving	EPA, VESA	EPA, VESA
MTBF	>50,000 hours Professional Protective glass	>50,000 hours Professional Protective glass