

## ALPHENDOSCOPES®

- Custom made “Sapphire objectives” guarantee lasting
- Autoclavability, scratch resistance and large field of view
- State of art ‘Rod Lens Systems’ for superior image quality
- High tech ‘Tubing Material’ more durable and scratch resistant and rust proof
- Optimized fiber optics provide more light



### LAPAROSCOPES FULL HD

Description	Part Number	Dia	Length	View Angle
Laparoscope	LPFH-05-310-0°	5mm	310	0°
Laparoscope	LPFH-05-310-30°	5mm	310	30°
Laparoscope	LPFH-10-340-0°	10mm	340	0°
Laparoscope	LPFH-10-340-30°	10mm	340	30°

### LAPAROSCOPES HD

Description	Part Number	Dia	Length	View Angle
Laparoscope	LPHD-10-315-0°	10mm	315	0°
Laparoscope	LPHD-10-315-30°	10mm	315	30°
Laparoscope Bariatric	LPHD-10-420-0°	10mm	420	0°
Laparoscope Bariatric	LPHD-10-420-30°	10mm	420	30°
Laparoscope Bariatric	LPHD-10-420-45°	10mm	420	45°

### LAPAROSCOPES STD

Description	Part Number	Dia	Length	View Angle
Laparoscope	LPST-05-300-0°	05mm	300	0°
Laparoscope	LPST-05-300-30°	05mm	300	30°
Laparoscope	LPST-65-300-0°	6.5mm	300	0°
Laparoscope	LPST-65-300-30°	6.5mm	300	30°
Laparoscope	LPST-10-315-0°	10mm	315	0°
Laparoscope	LPST-10-315-30°	10mm	315	30°

### SINOSCOPE / ARTHROSCOPE STD

Description	Part Number	Dia	Length	View Angle
Sinoscope	SINST-04-175-0°	04mm	175	0°
Sinoscope	SINST-04-175-30°	04mm	175	30°
Sinoscope	SINST-04-175-45°	04mm	175	45°
Sinoscope	SINST-04-175-70°	04mm	175	70°

### LAPAROSCOPE / CYSTOSCOPE / HYSTERSCOPE HD/ STD

Description	Part Number	Dia	Length	View Angle
Laparo/Cysto/Hystro	LCHST-29-300-0°	2.9mm	300	0°
Laparo/Cysto/Hystro	LCHST-29-300-30°	2.9mm	300	30°
Laparo/Cysto/Hystro	LCHHD-29-300-0°	2.9mm	300	0°
Laparo/Cysto/Hystro	LCHHD-29-300-30°	2.9mm	300	30°

### CYSTOSCOPE / HYSTERSCOPE HD/ STD

Description	Part Number	Dia	Length	View Angle
Cystoscope	CYSHD-04-302-0°	4mm	302	0°
Cystoscope	CYSHD-04-302-30°	4mm	302	30°
Cystoscope	CYSST-04-302-0°	4mm	302	0°
Cystoscope	CYSST-04-302-30°	4mm	302	30°

### LAPAROSCOPE SINGLE PUNCTURE [ offset ]

Description	Part Number	Dia	Length	View Angle
Laparoscope single puncture [6mm Dia Instrument channel]	LSP-10-270-0°	10mm	270	0°

\*Telescope Box  
TB-01/TB-02/TB-03



\*Protection Sheath  
TPS-10-31/TPS-10-42/TPS-4-18/TPS-45-31



TB-01 : Telescope box for Arthroscope / Sinoscopes  
TB-02 : Telescope box for Laparoscope  
TB-03 : Telescope box for Bariatric Scopes  
TB-04 : Telescope box for Cystoscope/Hystroscope

TPS-10-31 : Telescope Protection Sheath for 10mm Scopes  
TPS-10-42 : Telescope Protection Sheath for Bariatric Laparoscope  
TPS-4-18 : Telescope Protection Sheath for 4mm Sinoscope  
TPS-45-31 : Telescope Protection Sheath for 4/5mm Cystoscope & 5mm Laparoscope