


## ALPHKUT® EASY [AE-400A]

- Isolated RF Monopolar & Bipolar outputs.
- Well Suited for under water Tur-P & TU-BT.
- Maximum power output 400w.
- Separate footswitch for Monopolar & Bipolar.
- It is useful in open & Laparoscopic procedure.
- Portable & User Friendly design.
- Easy to mobility.
- Cost effective.



### DEVICE

Part No.	Description
AE-400A	Alphkut® Easy Analog [Electrosurgical Unit] With All Its Standard Accessories 400W (Monopolar & Bipolar)

### SUPPLIED ACCESSORIES

Part no.	Description
SSP-02	SS Patient plate
PPC-02	Patient plate cable
FSC-03-AE	Tri Pedal Foots witch
MCH-01	Monopolar chuck handle
MC-03	Mains Cord

### TECHNICAL DATA

MONOPOLAR CUT			
PURE CUT	400 Watt	500 OHM	1.5 CF
BLEND I	250 Watt	500 OHM	2.0 CF
BLEND II	200 Watt	500 OHM	2.5 CF
MONOPOLAR COAG			
SPRAY	120 Watt	500 OHM	7.0 CF
FORCE	150 Watt	500 OHM	6.0 CF
DESSICATE	100 Watt	500 OHM	5.0 CF
BIPOLAR			
BIPOLAR	80 Watt	100 OHM	1.5 CF
Power supply	170- 260 VAC		
Frequency	50 Hz		
Maximum Power	400W		
Classification	Protection Class 1, type Bf		
Dimension (W x H X D)	300 x 130 x 365 mm		
Weight	8kg		

SS Patient plate  
SSP-02



Patient plate cable  
PPC-02



Tri Pedal Footswitch  
FSC-03-AE



Mains cord  
MC-03



Monopolar chuck handle  
MCH-01



\*Bipolar cable  
BCS-04-03-VMF



\* Open bipolar tweezer  
OBT-07



\*Silicon Patient plate  
SPP-02

