


ALPHSORCE™ ENDO LED LIGHT [LSL-E300]

- Long life > 10000 hrs and more
- High Lumen output (Superior to 300w of Xenon)
- Instant on/off (no warm up or cool down)
- No change in colour temperature over entire color spectrum
- Does not emit ultraviolet and infrared light. Reduced heat transmission to the Endoscope
- LCD Display : Status and brightness control
- Overheat protection
- Consumes very less power
- No annual bulb change is necessary



DEVICE

Part No.	Description
LSL-E300	Alphsorce™ LED Light E300 For Endoscopy with Light cable socket,Olympus

SUPPLIED ACCESSORIES

Part No.	Description
MC-03	Mains cord

TECHNICAL DATA

Lamp	LED LSL-E300
Color Temperature	5600°k
Light guide connection	Olympus type , others on request
Insufflation Pump	Operate Air/water Function
Lumen Maintenance	90% After 10000hrs
Lumen Output	2700Lm
Ambient temperature	+10'C+/-40'C in operating mode
Cooling Method	Forced Air cooling side ventilation
Power Supply	100-240 AC
Frequency	50 to 60hz
Maximum Power	130W
Classification	Protection Class I, Type BF
Dimensions (W x H x D)	300 x 130 x 250 mm
Weight	7.5 Kg.

Mains cord
MC-03



* Air Bottle for olympus Endoscope
AB-01



*Fibre optic cable for Olympus LS
FOC-523-OS

