

**ALPHENDOTRAINER**

- Basic hand eye co-ordination and to learn minimal access surgery techniques (basic and advance)
- To learn intra and extra corporeal undo suturing and knot tying.
- Video Scopic Endo trainers are the most user friendly trainer which can be attached to any Monitor/TV.
- Built in illumination & high resolution camera
- 5 ports with silicon rubber gasket.

**DEVICE**

Part No.	Description
DSPT-1-C	Alphendo Trainer (Dome Shaped) with Camera, Light & Its Complete Accessories
DSPT-1-LK	Alphendo Trainer (Dome Shape) with Camera & Light Only
DSPT-1-L	Alphendo Trainer (Dome Shape) with Light
DSPT-1	Alphendo Trainer (Flat Shape) ( Box Only )
PTAS	Alphendo Trainer Accessories Set



**TECHNICAL DATA**

Power supply	100-240 V ac
Frequency	50-60Hz
Maximum Power	50W
Classification	Protection Class 1, type Bf
Dimension (W x H X D)	450mmX190mmX295mm
Weight	12kg


**SUPPLIED ACCESSORIES - DSPT-1-C**

Part no.	Description
APP-05-30-M3C	Plain Grasper 5mm -NDM
DMPT-05-30-M3C	Dissector Maryland 5mm - NDM
SCPT-05-30-M3C	Curve Scissor 5mm - NDM
NSPT-05-30	Needle Holder T.C. Axial Type straight 5mm - NDM
MC2-1	Mains 2 Pin Cable (for DSPT-1, DSPT-1-C)
RCA to RCA	Video cable
KPBR-05-40-12	Bhandarkar's knot Pusher (Ø=1.2mm) with Reducer

Bhandarkar's knot Pusher  
KPBR-05-40-12



Plain Grasper 5mm  
APP-05-30-M3C



Dissector Maryland 5mm  
DMPT-05-30-M3C



Curved Scissor 5mm  
SCPT-05-30-M3C



Needle Holder TC Axial Straight /  
Curve 5mm \*\*  
NSPT-05-30 /NCPT-05-30



Mains 2 Pins Cable  
MC2-1



Video cable  
RCA to RCA



\* Silicon Suture Pad  
SSP-4x4



\* Washer  
WET-05/WET-10



Further Accessories available on request