

ALPHSORCE™ HALOGEN LIGHT [LSH - 250]

- Special heat filters are used to eliminate thermal radiation.
- Electronically controlled light source with variable regulations, light output is very bright.
- Built in Voltage stabilizer and ventilation system adequate to dissipate heat from equipment.
- Instant Switch over to stbdy Bulb



DEVICE

Part No.	Description
LSH - 250	Alphsorce™ Halogen Cold Light - 250 W, Dual Bulb for Laparoscopy with Light cable socket, Storz

SUPPLIED ACCESSORIES

Part no.	Description
MC-03	Mains cord

TECHNICAL DATA

Lamp	Two optically centered halogen lamps with dichroic reflector PHL - 250w
Color Temperature	3200°k
Light guide connection	Karl Storz type , others on request
Life span of lamps	Min 50hrs at full power 100hrs at more or medium power.
Lumen Output	650 lumens
Ambient temperature	+10°C to +40°C in operating mode
Cooling Method	Forced Air cooling side ventilation
Power Supply	100-240 AC
Frequency	50 Hz
Maximum Power	400 VA
Classification	Protection Class I, Type BF
Dimensions (W x H x D)	226 x 120 x 270mm
Weight	6 Kg.

Mains cord
MC-03



* Halogen Lamps 24v/250w
PHL - 250w



*Fibre optic cable Armoured
FOC-523-SS



*Fibre optic cable Non-armoured
FXN-523-SS

