


ALPHIRRIGATION™ PUMP [IPex-1000]

- Intelligent leak compensation for optimal view an fluid balance
- Non contact pressure measurement
- Very easy to operate
- Permanent patient monitoring
- Quick start-up time.
- Easy to load tube set
- Precise control and regulation of pressure
- Small learning curve for surgeon and nurse



DEVICE

Part No.	Description
IPex -1000	Alphirrigation Pump [Roller Ball] [IPex-1000]

SUPPLIED ACCESSORIES

Part no.	Description
ITS-03	Irrigation Tube Set with IV Needle & Luer Lock
MC-03	Mains cord

TECHNICAL DATA

Insufflation medium	Saline / Glycerine
Pressure Range	5 to 400mmhg
Flow range	0.1 to 1LPM
Power Supply	170 to 240
Frequency	50 Hz
Maximum Power	150W
Classification	Protection Class I, Type BF
Dimensions (WxHxD)	300 x 320 x 138 mm
Weight	12 kgs.

Irrigation tubing set
ITS-03



Mains cord
MC-03

