


ALPHIRRIGATION™ PUMP [SIP-1600]

- Intelligent leak compensation for optimal view and fluid balance
- Non contact pressure measurement
- Good suction capability for surgery
- Very easy to operate
- Permanent patient monitoring
- Quick start-up time.
- Easy to load tube set
- Precise control and regulation of pressure & flow
- Precise control and regulation of Suction & flow
- Small learning curve for surgeon and nurse



DEVICE

Part No.	Description
SIP-1600	Alphirrigation Pump [Roller Ball] with Suction [SIP-1600]

SUPPLIED ACCESSORIES

Part no.	Description
ITS-03	Irrigation Tube Set with IV Needle & luer Lock
MC-03	Main cord
STS-03-06	Suction tubing set

Irrigation tubing set
ITS-03



Mains cord
MC-03



* Suction bottle
SB-4



Suction tubing set
STS-03-06



TECHNICAL DATA

Insufflation medium	Saline / Glycerine
Pressure Range	5 to 700 mmHg
Flow range	0.1 to 1.5lpm
Suction/ Vacuum	0 to - 560 mmHg
Power Supply	170 to 240 VAC
Frequency	50 Hz
Maximum Power	200w
Classification	Protection Class I, Type BF
Dimensions (WxHxD)	300 x 320 x 138 mm
Weight	12 kgs.