


**ALPHSORCE® LED LIGHT [LSL-300]**

- Long life > 10000 hrs and more
- High Lumen output ( Superior to 300W of Xenon)
- Instant on/off (no warm up or cool down)
- No change in color temperature over entire color spectrum
- Does not emit ultraviolet and infrared light. Reduced heat transmission to the Endoscope
- LCD Display shows Status and brightness control
- Overheat Protection and alarm
- Consumes very less power
- Fiber optic cable Sensing
- No annual bulb change is necessary



**DEVICE**

Part No.	Description
LSL-300	Alphsorce™ LED Light L300 for Laparoscopy with Light cable socket, Storz

**SUPPLIED ACCESSORIES**

Part no.	Description
MC-03	Mains cord

**TECHNICAL DATA**

Lamp	LED LSL-300
Color Temperature	5600°k
Light guide connection	Karl Storz type, others on request
Lumen Maintenance	90% After 10000hrs
Lumen Output	2700Lm
Ambient temperature	+10°C to +40°C in operating mode
Cooling Method	Forced Air cooling side ventilation
Power Supply	100-240 AC
Frequency	50 Hz
Maximum Power	100W
Classification	Protection Class I, Type BF
Dimensions (W x H x D)	300 x 130 x 250mm
Weight	5.5 Kgs.

Mains cord  
MC-03



\*Fibre optic cable Armoured  
FOC-523-SS



\*Fibre optic cable Non-armoured  
FXN-523-SS

