

TROLLEY

- Adjustables 5 shelves
- Eight point electrical board with socket & switch
- Co2 Cylinder cage
- Stainless Steel
- Lockable drawer
- Long service life
- Gliding wheels with locking brakes
- Stainless steel surface: easy cleaning, resistant to disinfectants agents

DEVICE

Part No.	Description
LP5-66-66-150	LAPAROSCOPIC TROLLEY



TECHNICAL DATA

Overall size	W 660mm x D 660mm x H 1498mm
shelves	5 Adjustables
Heavy duty castors	125 mm Dia